

FRANK OPTIC PRODUCTS®

.....
optische Technologien

*more than
just optics*

SideFire laser cable

for laser-induced therapy

FRANK OPTIC PRODUCTS® is a specialist and one of the leading manufacturers of medical probes and, developing and producing customised optical and fibre-optic components and systems for laser and minimal invasive surgery as well as in PDT.

As an OEM we supply **SideFire laser cables**, which are used as medical laser probes for the purpose of vaporisation, with customer-specific diverter angles at one end of the fibre and with customised laser connectors.

The **SideFire laser cables** are developed and supplied exclusively according to customer specifications and exclusively for their applications as medical laser probes for a variety of applications, laser types and laser outputs.

The **SideFire laser cables** made by FRANK OPTIC PRODUCTS® are distinguished itself by their great flexibility and high degree of stability.

The cable systems are available in a variety of diameters.

Standard diameter: $\varnothing_{\text{fibre core}}$:	200 μm to 1000 μm
Standard lengths:	Total length 1.5 m bis 5 m
Standard plug:	SMA-905
Holders:	customer-specific



All **SideFire laser cables** are supplied sterilised with CE and ETO in special biocompatible medical packaging, in conformity with EU Regulation 93/42 EEC.

FRANK OPTIC PRODUCTS GmbH optische Technologien, Heidelberger Straße 63-64, 12435 Berlin
phone +49(0)30530249-0 fax +49(0)30530249-21 e-mail info@fop-berlin.de www.frank-optic-products.de